

PATIENTS CAN  
**DO IT**  
WITHOUT  
FINGERSTICKS\*

## Introducing the FreeStyle Libre 14 day System

Glucose monitoring with a painless<sup>1</sup> scan instead of painful fingersticks\*

two components

### Innovative sensor technology

- Up to 14 days of accurate<sup>2</sup> continuous glucose readings day and night<sup>†</sup>
- Painless<sup>1</sup> sensor application
- No patient calibration or coding required
- Easy to wear<sup>1</sup> due to its small size
- Water resistant so you can swim<sup>‡</sup>, shower, or exercise
- Glucose readings one hour after application

### Compact, lightweight reader

- Shows the current glucose value, trend arrow, and 8 hours of glucose history
- Backlit color touch screen
- Stores 90 days of glucose data
- Compatible with FreeStyle Precision Neo test strips



Rx only

  
**FreeStyle**  
**Libre 14 day**  
FLASH GLUCOSE MONITORING SYSTEM

  
life. to the fullest.™

**Abbott**

\*Fingersticks are required for treatment decisions when you see Check Blood Glucose symbol, when symptoms do not match system readings, when you suspect readings may be inaccurate, or when you experience symptoms that may be due to high or low blood glucose.

†Based on the sensor being replaced once every 14 days, and scanned at least once every 8 hours.

‡Sensor is water-resistant in up to 1 meter (3 feet) of water. Do not immerse longer than 30 minutes.


REFERENCES: 1. Data on file. Abbott Diabetes Care. 2. FreeStyle Libre 14 day User's Manual. See back for Indications and Important Safety Information.

# FreeStyle Libre 14 day system basics:



## What is it?

FreeStyle Libre 14 day system is a **continuous glucose monitoring (CGM)** system consisting of a handheld reader and a sensor worn on the back of the upper arm



## How does it work?

The sensor has a thin, flexible filament **measuring glucose in interstitial fluid every minute**



## How do you use it?

Use the reader to perform a **painless<sup>1</sup> scan** instead of painful fingersticks\*

## How easy<sup>1</sup> is the sensor applied<sup>†</sup>?



**Apply the sensor** on the back of your upper arm



**Activate the sensor** with the reader to initiate the 1-hour sensor warm-up period



**Get a Reading** hold the reader within 1.5 inches of sensor to check glucose readings

## What does the reader show with each scan?

**A number** representing your current glucose reading

**A trend arrow** showing the direction your glucose is heading

**A trend graph** depicting the latest 8 hours



\*Fingersticks are required for treatment decisions when you see Check Blood Glucose symbol, when symptoms do not match system readings, when you suspect readings may be inaccurate, or when you experience symptoms that may be due to high or low blood glucose.

<sup>†</sup>For more detailed application instructions, please reference the User's Manual or Quick Start Guide. REFERENCE: 1. Data on file. Abbott Diabetes Care.

# What is included in the Reader and Sensor Kits?

## FreeStyle Libre 14 day Reader Kit

- 1 FreeStyle Libre 14 day Reader
- 1 USB Cable
- 1 Power Adapter
- User's Manual
- Quick Start Guide  
a guide to help get patients started quickly
- Treatment Decisions Guide  
a guide to help patients use reader information to make informed treatment decisions
- Interactive Tutorial USB Card  
a comprehensive electronic resource on how to operate the FreeStyle Libre 14 day system, including instructional videos



## FreeStyle Libre 14 day Sensor Kit

- 1 FreeStyle Libre 14 day Sensor Pack  
contains sensor to be worn for up to 14 days
- 1 Sensor Applicator  
applies the sensor to the back of the upper arm
- Alcohol Wipes
- Product Insert



The FreeStyle Libre 14 day system reader and sensors cannot be used interchangeably with the FreeStyle Libre system (10 day)

# FreeStyle Libre brand is the #1 CGM used worldwide<sup>1</sup>

A prescription is required to dispense the FreeStyle Libre 14 day system. There should be two parts to the prescription:

- 1 Annual Reader Prescription**  
NDC 57599-0002-00  
First time customers will need a reader with their first fill. FreeStyle Libre 14 day readers have a mean use life of three years.
- 2 Monthly Sensor Prescription**  
NDC 57599-0001-01  
The FreeStyle Libre 14 day sensor package contains 1 disposable sensor. A typical month should consist of 2 sensors (28 day supply)

- 1 RX: FreeStyle Libre 14 day Reader**  
**QTY: 1**
- 2 RX: FreeStyle Libre 14 day Sensor**  
**QTY: 2 Sensors a month**  
(28-day supply)  
**REFILLS: 12 (1-year supply)**

**Rx** PATIENT NAME \_\_\_\_\_  
ADDRESS \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
Prescription:  
  
SIGNATURE: \_\_\_\_\_  
DATE: \_\_\_\_\_

For illustrative purposes only.

**FreeStyle Libre 14 day Reader**  
NDC 57599-0002-00  
ABC 8: 10190381  
ABC 6: 782813  
Cardinal: 5479092  
McKesson: 2535201



**FreeStyle Libre 14 day Sensor**  
NDC 57599-0001-01  
ABC 8: 10190382  
ABC 6: 782817  
Cardinal: 5479084  
McKesson: 2535193

System components are sold in eaches  
Two (2) sensors for a 28-day supply

Learn more at [Pharmacy.FreeStyleLibre.us](https://Pharmacy.FreeStyleLibre.us)



**REFERENCE 1.** Data based on the number of users worldwide for FreeStyle Libre personal CGM compared to the number of users for other leading personal CGM brands.

#### Indications and Important Safety Information

The FreeStyle Libre 14 day Flash Glucose Monitoring System is a continuous glucose monitoring (CGM) device indicated for the management of diabetes in persons age 18 and older. It is designed to replace blood glucose testing for diabetes treatment decisions. The System detects trends and tracks patterns aiding in the detection of episodes of hyperglycemia and hypoglycemia, facilitating both acute and long-term therapy adjustments. Interpretation of the System readings should be based on the glucose trends and several sequential readings over time. The System is intended for single patient use and requires a prescription..

**CONTRAINDICATIONS:** Remove the sensor before MRI, CT scan, X-ray, or diathermy treatment.

**WARNINGS/LIMITATIONS:** Do not ignore symptoms that may be due to low or high blood glucose, hypoglycemic unawareness, or dehydration. Check sensor glucose readings with a blood glucose meter when Check Blood Glucose symbol appears, when symptoms do not match system readings, or when readings are suspected to be inaccurate. The FreeStyle Libre 14 day system does not have alarms unless the sensor is scanned, and the system contains small parts that may be dangerous if swallowed. The FreeStyle Libre 14 day system is not approved for pregnant women, persons on dialysis, or critically-ill population. Sensor placement is not approved for sites other than the back of the arm and standard precautions for transmission of blood borne pathogens should be taken. The built-in blood glucose meter is not for use on dehydrated, hypotensive, in shock, hyperglycemic-hyperosmolar state, with or without ketosis, neonates, critically-ill patients, or for diagnosis or screening of diabetes. Review all product information before use or contact Abbott toll-free 855-632-8658 or visit [www.FreeStyleLibre.us](https://www.FreeStyleLibre.us) for detailed indications for use and safety information. For full indications for use and safety information, visit <https://www.FreeStyleLibre.us/safety-information.html>.

FreeStyle, Libre, and related brand marks are trademarks of Abbott Diabetes Care Inc. in various jurisdictions. Other trademarks are property of their respective owners.  
©2018 Abbott. ADC-09306 v1.0 09/18